

Technical Information Sheet Product Number #1522 3M™ Double Coated Medical Tape

Effective: July 2013

Supersedes all previous specifications

Features and Benefits:

- Hypoallergenic pressure sensitive adhesive diecutting
- Heavy weight release liner for
- EtO & Gamma Sterilization Compatible
- Easy to Apply

Clear Film Carrier

- Adheres to a wide range of substrates

Composition:

Tape caliper, no liner	6.3 mil (0.16 mm) Double Coated Polyethylene Tape	
Carrier	3 mil (0.08 mm) Transparent Polyethylene Film	
Adhesive	Acrylate designed for medical/surgical use	
Release Liner	80 lb. Bleached Kraft paper, silicone release both sides (4.5 mil/0.11mm)	

Typical Properties⁺:

180° Peel Adhesion to Stainless Steel, Face [Open] side	28 oz/inch width (800 gm/25mm)
Liner Release	9 gm/inch width (9 gm/25 mm)

⁺ Properties listed based on very limited data and subject to change. Methods available upon request.

Roll Description: Tape supplied on 3 inch (76 mm) diameter cores.

Length	550 yards (503 m) ± 2% Maximum
Width	48 inch (122 cm) ± 1/32 in (1.0 mm)

Number of Splices: Splices on the tape are intended for continuous processing.

Roll Width	Maximum Number of Splices per Roll
≤ 8 inch (20 cm)	2 per 100 yd (91 m)
	3 per 200 yd (183 m)
	6 per 550 yd (503 m)

Information on the splice frequency in roll widths greater than 8 inch (20 cm) is available upon request.

Packaging and Recommended Storage

Product as supplied in original packaging will maintain stated properties for a period of three (3) years from date stamped on shipping container when stored at temperatures of 50-80°F (10-27°C) and a relative humidity of 40-60 percent.

Please see reverse side for important product, safety and warranty information.

Product and Safety Information: User is solely responsible for determining the suitability of 3M samples and products for the intended use including any necessary safety or toxicity assessment. 3M will provide Material Safety Data Sheets and results of toxicity testing upon request. In every case before using any product in full scale production users should conduct their own tests to determine to their own satisfaction whether the product is of acceptable quality and is suitable for their particular purposes under their own operating conditions.

Notice: Nothing contained herein shall be construed to imply the nonexistence of any relevant patents or to constitute a permission, inducement or recommendations to practice any invention covered by any patent, without authority from the owners of this patent.

All statements, technical information and recommendations herein are based on tests 3M believes to be reliable, but the accuracy or completeness thereof is not guaranteed. 3M warrants only that products will meet 3M's specifications at the time of shipment to the customer. 3M does not offer any other warranty and does not warrant the performance, safety or such other characteristics of Products in combination with other materials. 3M specifically DOES NOT warrant Products for any intended or unintended uses (whether or not foreseeable); for compatibility or suitability with other components or compatibility with any methods of manufacture or conversion. The foregoing warranty is made in lieu of all other warranties, expressed or implied, including the implied warranties of merchantability, fitness for a particular purpose and freedom from non-infringement.

Limitation of Remedies: If products are proven not to meet 3M's specifications, the sole and exclusive remedy available and 3M's only obligation shall be, at 3M's option, to replace such quantity of Products which are proven out of specification or to refund the purchase price paid for Products.

Limitations of Liabilities: The remedies provided herein are exclusive remedies against 3M for any alleged or actual nonconformance to specifications or defect or other failure in products or for 3M's performance of its supply obligations. Under no circumstances is 3M liable for any direct, indirect, incidental, special or consequential damages (including lost profits) in any way related to the product under any theory of law including, but not limited to, negligence and strict liability.

Ordering Information

Call us at Customer Service to place an order: 800-742-1994 (U.S.). Visit our website: www.3M.com/medicalspecialties for product and services information, news about conferences we will be attending, new product highlights or to make a direct inquiry.

To have a sales representative contact you, to request samples or clinical and safety summaries, please contact us at 3M HELPLINE 800-228-3957 (U.S.) or for international inquiries, please contact your local country representative. Our 3M Medical Specialties subsidiary contacts are listed below.

International Locations

Argentina **Dominican Republic** Kazakhstan **Poland** Thailand 3M Argentina S.A.C.I.F.I.A. 3M Poland Sp.z o.o. 3M Dominicana 3M Representation Office 3M Thailand Ltd. Tel.: 54-11-4339-2400 Tel: (809) 530-6560 X499 Tel: +7 3272 509 944 Tel: (48 22) 739 60 00 Tel: (66)2 260 8577 Ext.175 Fax: 54-11-4339-2640 Fax: (809) 530-2960 Fax: +7 3272 509 573 Fax: (48 22) 739 60 01 Fax: (66)2 261 7535 Australia Finland Portugal Taiwan Korea 3M Health Care Suomen 3M Oy 3M Korea Ltd. 3M Portugal 3M Taiwan Ltd. Tel: 1 300 363 878 Tel.: 358-9-52-521 Tel: 82-2-3771-4286 Tel: 351-21-3134500 Tel: 886-2-2704-9011 Fax: 1 800 060 888 Fax: 358-9-512-29-44 Fax: 82-2-786-2825 Fax: 351-21-3134680 Fax: 886-2-2706-0355 **Belgium** France Malaysia Russia Turkey 3M Belgium N.V./S.A. Laboratoires 3M Santé 3M Malaysia Sdn. Bhd. 3M Russia 3M Turkey Tel: 6-03-7806 2888 Tel: 9 (0212) 3507777 Tel: 32-2-722-5111 Tel: 33-1-3031-8376 Tel: 7 (495) 784-74-74 Fax: 32-2-720-0225 Fax: 33-1-3031-8378 Fax: 6-03-7806 2902 Fax: 7 (495) 784-74-75 Fax: 9 (0212) 2821742 Brazil Germany Singapore Ukraine 3M do Brasil Ltda. 3M México, S.A. de C.V. 3M Technologies (S) Pte. Ltd. 3M Ukraine 3M Medica Tel: 55-19-3838-6661 Tel: 49 2131 14 40 00 Tel: +52-55-5270-0400 Tel: 65-6450-8866 Fax: 55-19-3838-7449 Fax: 49 2131 14 49 99 Fax: +52-55-5270-0433 Fax: 65-6455-2130 Canada Hong Kong Middle East Sri Lanka 3M Canada 3M Hong Kong Ltd. 3M Gulf Ltd. 3M Lanka (Pvt.) Ltd. Tel: 971-4-3670-777 Tel: (94) 11 2785701

Tel: 800 364-3577 Tel: (852) 2806 6111 Fax: (852) 2234 6044 Fax: 800 341-4630

Chile Hungary 3M Chile S.A. 3M Hungária Kft. Tel: (56-2) 4103000 Tel: (36-1) 270-7777 Fax: (56-2) 4103400 Fax: (36-1) 320-0951

China India 3M China Ltd. 3M India Ltd. Tel: 86-21-6275-3535 Tel: +91-80-2223 1414 Fax: 86-21-5208-2205 Fax: +91-80-22231431

Colombia Indonesia 3M Colombia S.A. PT 3M Indonesia Tel: 57-1 4101666 Tel: (62-21) 520 3401 Fax: 57-1 4161677 Fax: (62-21) 520 3106

Czech Republic 3M Cesko s.r.o. 3M Italia SPA Tel: +420-261 380 111 Tel: +39 2 7035-2406 Fax: +420-261 380 110 Fax: +39 2 7035-2484

Denmark 3M Health Care Ltd. 3M a/sTel: 81-3-3709-8289 Tel: 45-43-48-0100 Fax: 45-43-96-8596 Fax: 81-3-3709-8754

New Zealand 3M New Zealand Tel: 64-9-444 4760 Fax: 64-9-443 7885

Fax: 971-4-3670-699

The Netherlands

3M Nederland B.V.

Tel: 31-715-450-450

Fax: 31-715-450-212

Pern 3M Peru S.A. Tel: 511-224-2728 Fax: 511-224-3171

Philippines 3M Philippines Tel: (632) 8133781 to 95 Fax: (632) 8145873 to 74 Tel: +380 44 490 57 77 Fax: +380 44 490 57 75 **United Kingdom & Ireland** 3M Health Care Ltd. Tel: 44-1509-611-611

Fax: 44-1509-237-288

Spain Venezuela 3M España, S.A. 3M Manufacturera Venezuela Tel: 34-91-321-6000 Tel: 58-241-8391911 Fax: 34-91-321-6002 Fax: 58-241-8391928

Viet Nam Sweden 3M Svenska AB 3M Viet Nam Ltd. Tel: +46-8-92 21 00 Tel: 84 8 640986 Fax: +46-8-92 24 22 Fax: 84 8 641024

Switzerland & Eastern Europe

3M (Schweiz) AG Tel: +41 44 724 92 44 Fax: +41 44 724 94 80

Fax: (94) 11 4410083-4



3M Medical Specialties 3M Center, Building 275-5W-05 St. Paul. MN 55144-1000 USA Tel.: 800 228-3957

