

Technical Information Sheet Product Number 1523 3MTM Blush Tan Polyethylene Single Coated Medical Tape

Effective: July 2013

Supersedes all previous specifications

Features and Benefits:

- Hypoallergenic Pressure Sensitive Adhesive

ETO & Gamma Sterilization CompatibleGood Adhesion to Skin

- Conformable

- Opaque Blush Tan

Easy to Apply

- Provides Barrier from External Contamination - No Na

- No Natural Latex Rubber

Composition:

Tape Caliper	5 mil (0.13 mm) Tan Polyethylene Film Tape
Backing	3.6 mil (0.09 mm) Blush Tan Polyethylene
Adhesive	Tackified Acrylic, designed for medical/surgical applications
Release Liner	63 lb. Poly-coated Kraft paper, silicone on one side (4.9 mils/0.12 mm)

Typical Properties⁺:

Adhesion to Steel, 180° Peel	41 oz/inch width (11.3 N/25 mm width) (1.2 Kg/25 mm width)
Liner Release	20 gm/inch width (20 gm/25 mm width)

⁺³M test methods are available upon request.

Roll Description: Tape supplied on 3 inch (76 mm) diameter cores.

Length, Maximum	500 yards (457 m) ± 2%
Width, Maximum	48 inch (122 cm) ± 1/32 in (1.0 mm)

Splices: Splices on the tape are intended for continuous processing.

Roll Width	Maximum Number of Splices per Roll
≤ 8 inch (20 cm)	2 per 100 yd (91 m)
	3 per 200 yd (183 m)
	8 per 500 yd (457 m)

Information on the splice frequency in roll widths greater than 8 inch (20 cm) is available upon request.

Packaging and Recommended Storage:

Product as supplied in original packaging will maintain stated properties for a period of two years from date stamped on shipping container when stored at temperatures of 50-80°F (10-27°C) and a relative humidity between 40-60 percent.

Please see reverse side for important product, safety and warranty information.

Product and Safety Information: User is solely responsible for determining the suitability of 3M samples and products for the intended use including any necessary safety or toxicity assessment. 3M will provide Material Safety Data Sheets and results of toxicity testing upon request. In every case before using any product in full scale production users should conduct their own tests to determine to their own satisfaction whether the product is of acceptable quality and is suitable for their particular purposes under their own operating conditions.

Notice: Nothing contained herein shall be construed to imply the nonexistence of any relevant patents or to constitute a permission, inducement or recommendations to practice any invention covered by any patent, without authority from the owners of this patent.

All statements, technical information and recommendations herein are based on tests 3M believes to be reliable, but the accuracy or completeness thereof is not guaranteed. 3M warrants only that products will meet 3M's specifications at the time of shipment to the customer. 3M does not offer any other warranty and does not warrant the performance, safety or such other characteristics of Products in combination with other materials. 3M specifically DOES NOT warrant Products for any intended or unintended uses (whether or not foreseeable); for compatibility or suitability with other components or compatibility with any methods of manufacture or conversion. The foregoing warranty is made in lieu of all other warranties, expressed or implied, including the implied warranties of merchantability, fitness for a particular purpose and freedom from non-infringement.

Limitation of Remedies: If products are proven not to meet 3M's specifications, the sole and exclusive remedy available and 3M's only obligation shall be, at 3M's option, to replace such quantity of Products which are proven out of specification or to refund the purchase price paid for Products.

Limitations of Liabilities: The remedies provided herein are exclusive remedies against 3M for any alleged or actual nonconformance to specifications or defect or other failure in products or for 3M's performance of its supply obligations. Under no circumstances is 3M liable for any direct, indirect, incidental, special or consequential damages (including lost profits) in any way related to the product under any theory of law including, but not limited to, negligence and strict liability.

Ordering Information

Call us at Customer Service to place an order: 800-742-1994 (U.S.). Visit our website: www.3M.com/medicalspecialties for product and services information, news about conferences we will be attending, new product highlights or to make a direct inquiry.

To have a sales representative contact you, to request samples or clinical and safety summaries, please contact us at 3M HELPLINE 800-228-3957 (U.S.) or for international inquiries, please contact your local country representative. Our 3M Medical Specialties subsidiary contacts are listed below.

International Locations

Argentina 3M Argentina S.A.C.I.F.I.A. Tel.: 54-11-4339-2400 Fax: 54-11-4339-2640
Australia 3M Health Care Tel: 1 300 363 878
Fax: 1 800 060 888
Belgium 3M Belgium N.V./S.A. Tel: 32-2-722-5111 Fax: 32-2-720-0225

3M do Brasil Ltda.
Tel: 55-19-3838-6661
Fax: 55-19-3838-7449
Canada

Brazil

3M Canada Tel: 800 364-3577

China

Chile
3M Chile S.A.
Tel: (56-2) 4103000
Fax: (56-2) 4103400

Fax: 800 341-4630

· · · · · · · · · · · · · · · · · · ·
3M China Ltd.
Tel: 86-21-6275-3535
Fax: 86-21-5208-2205

Colombia
3M Colombia S.A.
Tel: 57-1 4101666
Fax: 57-1 4161677

Czech Republic

Cheen republic
3M Cesko s.r.o.
Tel: +420-261 380 111
Fax: +420-261 380 110

Denmark	
3M a/s	
Tel: 45-43-48-0100	
Fav: 15-13-96-8596	

Dominican Republic

3M Dominicana
Tel: (809) 530-6560 X499
Eart (900) 520 2060

Finland
Suomen 3M Oy
Tel.: 358-9-52-521
Fax: 358-9-512-29-44

France
Laboratoires 3M Santé
Tel: 33-1-3031-8376
Fax: 33-1-3031-8378

Germany
3M Medica
Tel: 49 2131 14 40 00
Fax: 49 2131 14 49 99

Hong Kong
3M Hong Kong Ltd.
Tel: (852) 2806 6111
Eart (952) 2224 604/

Hungary
3M Hungária Kft.
Tel: (36-1) 270-7777
Fax: (36-1) 320-0951

Indonesia
PT 3M Indonesia
Tel: (62-21) 520 3401
Fax: (62-21) 520 3106

Israel
3M Israel Ltd.
Tel: 972-9-9615000
Fax: 972-9-9615050

Italy
3M Italia SPA
Tel: +39 2 7035-2406
Fax: +39 2 7035-2484

Japan
3M Health Care Ltd.
Tel: 81-3-3709-8289
Fax: 81-3-3709-8754

Kazakhstan
3M Representation Office
Tel: +7 3272 509 944
Fax: +7 3272 509 573

Korea
3M Korea Ltd.
Tel: 82-2-3771-4286
Fax: 82-2-786-2825

Malaysia
3M Malaysia Sdn. Bhd.
Tel: 6-03-7806 2888
Fax: 6-03-7806 2902

Mexico
3M México, S.A. de C.V
Tel: +52-55-5270-0400
Fax: +52-55-5270-0433

Middle East
3M Gulf Ltd.
Tel: 971-4-3670-777
Fax: 971-4-3670-699

The Netherlands
3M Nederland B.V.
Tel: 31-715-450-450
Fax: 31-715-450-212

New Zealand
3M New Zealand
Tel: 64-9-477 4040
Fax: 64-9-477 6692

Peru
3M Peru S.A.
Tel: 511-224-2728
Fax: 511-224-3171

Philippines
3M Philippines
Tel: (632) 8133781 to 95
Fax: (632) 8145873 to 74

Poland
3M Poland Sp.z o.o.
Tel: (48 22) 739 60 00
Fax: (48 22) 739 60 01

Portugal
3M Portugal
Tel: 351-21-3134500
Fax: 351-21-3134680

Russia	
3M Russia	
Tel: 7 (495) 784-74	-74
Fax: 7 (495) 784-74	I-75

Singapore	
3M Technologies (S) P	te. Ltd
Tel: 65-6450-8866	
Fax: 65-6455-2130	

Sri Lanka
3M Lanka (Pvt.) Ltd.
Tel: (94) 11 2785701

Fax: (94) 11 4410083-4

Spain
3M España, S.A.
Tel: 34-91-321-6000
Fax: 34-91-321-6002

Sweden	
3M Svenska AB	
Tel: +46-8-92 21 00	

3M Svenska AB
Tel: +46-8-92 21 00
Fax: +46-8-92 24 22

Switzerland & Eastern Europe

3M (Schweiz) AG
Tel: +41 44 724 92 44
Fax: +41 44 724 94 80



Thailand 3M Thailand Ltd. Tel: (66)2 260 8577 Ext.175 Fax: (66)2 261 7535

Taiwan 3M Taiwan Ltd. Tel: 886-2-2704-9011 Fax: 886-2-2706-0355

Turkey 3M Turkey

Ukraine 3M Ukraine

Tel: 9 (0212) 3507777 Fax: 9 (0212) 2821742

Tel: +380 44 490 57 77 Fax: +380 44 490 57 75

Tel: 44-1509-611-611

Fax: 44-1509-237-288

United Kingdom & Ireland 3M Health Care Ltd.

Fax: 58-241-8391928 Viet Nam 3M Viet Nam Ltd. Tel: 84 8 640986

Fax: 84 8 641024



3M Medical Specialties

3M Center, Building 275-5W-05 St. Paul. MN 55144-1000 USA

Tel.: 800 228-3957

