

# Technical Information Sheet Product Number 1530L 3M™ Medical White Rayon Nonwoven Tape on Liner

Effective: August 2013

Supersedes all previous specifications

#### **Features and Benefits:**

- Hypoallergenic Pressure Sensitive Adhesive - ETO and Gamma Sterilization Compatible

Easy to Tear by Hand - High MVTR

Breathable - Very Gentle, No Natural Rubber Latex

#### **Composition:**

Tape Caliper	5.5 mil (0.14 mm) White Nonwoven Rayon Tape	
Adhesive	Non-tackified Acrylate, designed for medical/surgical applications	
Release Liner	63 lb. Poly-coated Kraft paper, silicone on one side (4.6 mils/0.12 mm)	

# **Typical Properties**<sup>+</sup>:

Adhesion to Stainless Steel	17 ounce/inch width (480 gms/25 mm) (4.7 N/25 mm)		
Tensile Strength, at break, min.	5 lbs/inch width (2.26 Kg/25mm width) (22.2 N/25 mm)		
Moisture Vapor Transmission	$4200 \text{ gms/m}^2/24 \text{ hrs}$		
Removal from Liner	34 gm/inch (34 gm/25 mm)		

<sup>+3</sup>M test methods are available upon request.

# **Roll Description:** Tape supplied on 3 inch (76 mm) diameter cores.

Length	650 yards (594 meters) ± 2% Maximum	
Width	48 inch (122 cm) ± 1/16 in (1.5 mm) Maximum	

## **Splices:** Splices on the tape are intended for continuous processing.

Roll Width	Maximum Number of Splices per Roll	
≤ 8 inch (20 cm)	3 per 100 yd (91 m)	
	8 per 300 yd (274 m)	
	11 per 650 yd (594 m)	

Information on the splice frequency in roll widths greater than 8 inch (20 cm) is available upon request.

## **Packaging and Recommended Storage**

Product as supplied in original packaging will maintain stated properties for a period of three [3] years from date stamped on shipping container when stored at temperatures of 50-80°F (10-27°C) and a relative humidity of 40-60 percent.

Please see reverse side for important product, safety and warranty information.

Product and Safety Information: User is solely responsible for determining the suitability of 3M samples and products for the intended use including any necessary safety or toxicity assessment. 3M will provide Material Safety Data Sheets and results of toxicity testing upon request. In every case before using any product in full scale production users should conduct their own tests to determine to their own satisfaction whether the product is of acceptable quality and is suitable for their particular purposes under their own operating conditions.

Notice: Nothing contained herein shall be construed to imply the nonexistence of any relevant patents or to constitute a permission, inducement or recommendations to practice any invention covered by any patent, without authority from the owners of this patent.

All statements, technical information and recommendations herein are based on tests 3M believes to be reliable, but the accuracy or completeness thereof is not guaranteed. 3M warrants only that products will meet 3M's specifications at the time of shipment to the customer. 3M does not offer any other warranty and does not warrant the performance, safety or such other characteristics of Products in combination with other materials. 3M specifically DOES NOT warrant Products for any intended or unintended uses (whether or not foreseeable); for compatibility or suitability with other components or compatibility with any methods of manufacture or conversion. The foregoing warranty is made in lieu of all other warranties, expressed or implied, including the implied warranties of merchantability, fitness for a particular purpose and freedom from non-infringement.

Limitation of Remedies: If products are proven not to meet 3M's specifications, the sole and exclusive remedy available and 3M's only obligation shall be, at 3M's option, to replace such quantity of Products which are proven out of specification or to refund the purchase price paid for Products.

Limitations of Liabilities: The remedies provided herein are exclusive remedies against 3M for any alleged or actual nonconformance to specifications or defect or other failure in products or for 3M's performance of its supply obligations. Under no circumstances is 3M liable for any direct, incidental, special or consequential damages (including lost profits) in any way related to the product under any theory of law including, but not limited to, negligence and strict

# **Ordering Information**

Call us at Customer Service to place an order: 800-742-1994 (U.S.). Visit our website: www.3M.com/medicalspecialties for product and services information, news about conferences we will be attending, new product highlights or to make a direct inquiry.

To have a sales representative contact you, to request samples or clinical and safety summaries, please contact us at 3M HELPLINE 800-228-3957 (U.S.) or for international inquiries, please contact your local country representative. Our 3M Medical Specialties subsidiary contacts are listed below.

International Locations						
	Argentina 3M Argentina S.A.C.I.F.I.A. Tel.: 54-11-4339-2400 Fax: 54-11-4339-2640	<b>Dominican Republic</b> 3M Dominicana Tel: (809) 530-6560 X499 Fax: (809) 530-2960	<b>Japan</b> 3M Health Care Ltd. Tel: 81-3-3709-8289 Fax: 81-3-3709-8754	Poland 3M Poland Sp.z o.o. Tel: (48 22) 739 60 00 Fax: (48 22) 739 60 01	<b>Thailand</b> 3M Thailand Ltd. Tel: (66)2 260 8577 Ext.175 Fax: (66)2 261 7535	
	Australia 3M Health Care Tel: 1 300 363 878 Fax: 1 800 060 888	Finland Suomen 3M Oy Tel.: 358-9-52-521 Fax: 358-9-512-29-44	<b>Kazakhstan</b> 3M Representation Office Tel: +7 3272 509 944 Fax: +7 3272 509 573	Portugal 3M Portugal Tel: 351-21-3134500 Fax: 351-21-3134680	<b>Taiwan</b> 3M Taiwan Ltd. Tel: 886-2-2704-9011 Fax: 886-2-2706-0355	
	<b>Belgium</b> 3M Belgium N.V./S.A. Tel: 32-2-722-5111 Fax: 32-2-720-0225	France Laboratoires 3M Santé Tel: 33-1-3031-8376 Fax: 33-1-3031-8378	Korea 3M Korea Ltd. Tel: 82-2-3771-4286 Fax: 82-2-786-2825	<b>Russia</b> 3M Russia Tel: 7 (495) 784-74-74 Fax: 7 (495) 784-74-75	<b>Turkey</b> 3M Turkey Tel: 9 (0212) 3507777 Fax: 9 (0212) 2821742	
	Brazil 3M do Brasil Ltda. Tel: 55-19-3838-6661 Fax: 55-19-3838-7449	<b>Germany</b> 3M Medica Tel: 49 2131 14 40 00 Fax: 49 2131 14 49 99	Malaysia 3M Malaysia Sdn. Bhd. Tel: 6-03-7806 2888 Fax: 6-03-7806 2902	Singapore 3M Technologies (S) Pte. Ltd. Tel: 65-6450-8866 Fax: 65-6455-2130	Ukraine 3M Ukraine Tel: +380 44 490 57 77 Fax: +380 44 490 57 75	
	<b>Canada</b> 3M Canada Tel: 800 364-3577 Fax: 800 341-4630	Hong Kong 3M Hong Kong Ltd. Tel: (852) 2806 6111 Fax: (852) 2234 6044	<b>Mexico</b> 3M México, S.A. de C.V. Tel: +52-55-5270-0400 Fax: +52-55-5270-0433	<b>Sri Lanka</b> 3M Lanka (Pvt.) Ltd. Tel: (94) 11 2785701 Fax: (94) 11 4410083-4	United Kingdom & Ireland 3M Health Care Ltd. Tel: 44-1509-611-611 Fax: 44-1509-237-288	
	Chile	Hungary	Middle East	Spain	Venezuela	

3M Gulf Ltd.

Tel: 971-4-3670-777

Fax: 971-4-3670-699

The Netherlands

New Zealand

3M Nederland B.V.

Tel: 31-715-450-450

Fax: 31-715-450-212

Chile Hungary 3M Chile S.A. 3M Hungária Kft. Tel: (56-2) 4103000 Tel: (36-1) 270-7777 Fax: (56-2) 4103400 Fax: (36-1) 320-0951

China India 3M China Ltd. 3M India Ltd. Tel: 86-21-6275-3535 Tel: +91-80-2223 1414 Fax: 86-21-5208-2205 Fax: +91-80-22231431

Colombia Indonesia 3M Colombia S.A. PT 3M Indonesia Tel: 57-1 4101666 Tel: (62-21) 520 3401 Fax: 57-1 4161677 Fax: (62-21) 520 3106

Czech Republic 3M Cesko s.r.o. 3M Israel Ltd. Tel: +420-261 380 111 Tel: 972-9-9615000 Fax: +420-261 380 110 Fax: 972-9-9615050

Denmark 3M a/s3M Italia SPA Tel: 45-43-48-0100 Tel: +39 2 7035-2406 Fax: 45-43-96-8596 Fax: +39 2 7035-2484

3M New Zealand Tel: 64-9-477 4040 Fax: 64-9-477 6692 3M Peru S.A. Tel: 511-224-2728

Fax: 511-224-3171 **Philippines** 3M Philippines Tel: (632) 8133781 to 95 Fax: (632) 8145873 to 74 Switzerland & Eastern Europe 3M (Schweiz) AG Tel: +41 44 724 92 44 Fax: +41 44 724 94 80



3M Manufacturera Venezuela

Tel: 58-241-8391911

Fax: 58-241-8391928

3M Viet Nam Ltd.

Tel: 84 8 640986

Fax: 84 8 641024

Viet Nam

3M España, S.A.

3M Svenska AB

Tel: +46-8-92 21 00

Fax: +46-8-92 24 22

Sweden

Tel: 34-91-321-6000

Fax: 34-91-321-6002

**3M Medical Specialties** 3M Center, Building 275-5W-05 St. Paul. MN 55144-1000 USA Tel.: 800 228-3957