

Preliminary Technical Information Sheet

3MTM Medical Specialty Tape

1597 Two-in-One DC Polyester Tape (formerly MSX-7018NP)

Effective Date: April 2015

Supersedes all previous information sheets

Product Description:

#1597 is differential double-coated polyester film tape that uses synthetic rubber and tackified acrylic adhesives that were designed for medical applications. The tape is supplied on a double-coated 60# Bleached Kraft paper release liner covering the acrylic adhesive side of the tape.

Features and Benefits:

- Two Adhesive Systems
- Tackified acrylic adhesive for skin adhesion
- Synthetic rubber adhesive for low surface energy (LSE) plastics
- No Natural Rubber Latex
- Medical and Retail Device Applications
- EtO and Gamma compatible

Prototype Composition:

Tape caliper, no liner	3.9 mils (0.1 mm) Synthetic rubber/Acrylic '2-in-One' DC Tape		
Carrier	0.92 mils (0.02 mm) Transparent Polyester Film		
Adhesives, Open (Face) Side	Synthetic Rubber-based, designed for medical use		
Paper (Back) Side	Tackified Acrylic, developed for medical and surgical use		
Dalaaga I in an	60 lb. Printed Bleached Kraft paper, silicone release		
Release Liner	(3.5mil/0.09mm)		

Prototype Properties¹:

Adhesion to Stainless Steel, SR Side	175 oz/inch width (5 kg/25.4 mm)(48.6 Newtons)
Acr Side	40 oz/inch width (1.1 kg/25.4 mm)(11.1 Newtons)
Liner Release, Preliminary (Acr side)	10 gm/inch width (10 gm/25 mm)

¹Properties listed based on very limited data and subject to change.

Prototype Roll Description: Tape supplied on 3-inch (76 mm) diameter cores.

Length	700 yards (640 m) ± 2% Maximum
Width	48 inch (122 cm) \pm 1/32 in (1.0 mm)

Prototype Shelf Life and Storage Conditions:

Product as supplied in original packaging will maintain stated test properties for a period of up to two years when stored at temperatures of 50-80°F (10-27°C) and a relative humidity between 40-60 percent.

Please see reverse side for important product, safety and warranty information.

Product and Safety Information: User is solely responsible for determining the suitability of 3M samples and products for the intended use including any necessary safety or toxicity assessment. 3M will provide Material Safety Data Sheets and results of toxicity testing upon request. In every case before using any product in full scale production users should conduct their own tests to determine to their own satisfaction whether the product is of acceptable quality and is suitable for their particular purposes under their own operating conditions.

Notice: Nothing contained herein shall be construed to imply the nonexistence of any relevant patents or to constitute a permission, inducement or recommendations to practice any invention covered by any patent, without authority from the owners of this patent.

All statements, technical information and recommendations herein are based on tests 3M believes to be reliable, but the accuracy or completeness thereof is not guaranteed. 3M warrants only that products will meet 3M's specifications at the time of shipment to the customer. 3M does not offer any other warranty and does not warrant the performance, safety or such other characteristics of Products in combination with other materials. 3M specifically DOES NOT warrant Products for any intended or unintended uses (whether or not foreseeable); for compatibility or suitability with other components or compatibility with any methods of manufacture or conversion. The foregoing warranty is made in lieu of all other warranties, expressed or implied, including the implied warranties of merchantability, fitness for a particular purpose and freedom from non-infringement.

Limitation of Remedies: If products are proven not to meet 3M's specifications, the sole and exclusive remedy available and 3M's only obligation shall be, at 3M's option, to replace such quantity of Products which are proven out of specification or to refund the purchase price paid for Products.

Limitations of Liabilities: The remedies provided herein are exclusive remedies against 3M for any alleged or actual nonconformance to specifications or defect or other failure in products or for 3M's performance of its supply obligations. Under no circumstances is 3M liable for any direct, indirect, incidental, special or consequential damages (including lost profits) in any way related to the product under any theory of law including, but not limited to, negligence and strict liability.

Ordering Information

Call us at Customer Service to place an order: 800-742-1994 (U.S.). Visit our website: www.3M.com/medicalspecialties for product and services information, news about conferences we will be attending, new product highlights or to make a direct inquiry.

To have a sales representative contact you, to request samples or clinical and safety summaries, please contact us at 3M HELPLINE 800-228-3957 (U.S.) or for international inquiries, please contact your local country representative. Our 3M Medical Specialties subsidiary contacts are listed below.

h

Czech Republic 3M Cesko s.r.o.

Denmark

3M a/s

Tel: +420-261 380 111

Fax: +420-261 380 110

Tel: 45-43-48-0100

Fax: 45-43-96-8596

3M Israel Ltd.

3M Italia SPA

Tel: 972-9-9615000

Fax: 972-9-9615050

Tel: +39 2 7035-2406

Fax: +39 2 7035-2484

International Locations							
Argentina 3M Argentina S.A.C.I.F.I.A. Tel.: 54-11-4339-2400 Fax: 54-11-4339-2640	Dominican Republic 3M Dominicana Tel: (809) 530-6560 X499 Fax: (809) 530-2960	Japan 3M Health Care Ltd. Tel: 81-3-3709-8289 Fax: 81-3-3709-8754	Poland 3M Poland Sp.z o.o. Tel: (48 22) 739 60 00 Fax: (48 22) 739 60 01	Thailand 3M Thailand Ltd. Tel: (66)2 260 8577 Ext.175 Fax: (66)2 261 7535			
Australia 3M Health Care Tel: 1 300 363 878 Fax: 1 800 060 888	Finland Suomen 3M Oy Tel.: 358-9-52-521 Fax: 358-9-512-29-44	Kazakhstan 3M Representation Office Tel: +7 3272 509 944 Fax: +7 3272 509 573	Portugal 3M Portugal Tel: 351-21-3134500 Fax: 351-21-3134680	Taiwan 3M Taiwan Ltd. Tel: 886-2-2704-9011 Fax: 886-2-2706-0355			
Belgium 3M Belgium N.V./S.A. Tel: 32-2-722-5111 Fax: 32-2-720-0225	France Laboratoires 3M Santé Tel: 33-1-3031-8376 Fax: 33-1-3031-8378	Korea 3M Korea Ltd. Tel: 82-2-3771-4286 Fax: 82-2-786-2825	Russia 3M Russia Tel: 7 (495) 784-74-74 Fax: 7 (495) 784-74-75	Turkey 3M Turkey Tel: 9 (0212) 3507777 Fax: 9 (0212) 2821742			
Brazil 3M do Brasil Ltda. Tel: 55-19-3838-6661 Fax: 55-19-3838-7449	Germany 3M Medica Tel: 49 2131 14 40 00 Fax: 49 2131 14 49 99	Malaysia 3M Malaysia Sdn. Bhd. Tel: 6-03-7806 2888 Fax: 6-03-7806 2902	Singapore 3M Technologies (S) Pte. Ltd. Tel: 65-6450-8866 Fax: 65-6455-2130	Ukraine 3M Ukraine Tel: +380 44 490 57 77 Fax: +380 44 490 57 75			
Canada 3M Canada Tel: 800 364-3577 Fax: 800 341-4630	Hong Kong 3M Hong Kong Ltd. Tel: (852) 2806 6111 Fax: (852) 2234 6044	Mexico 3M México, S.A. de C.V. Tel: +52-55-5270-0400 Fax: +52-55-5270-0433	Sri Lanka 3M Lanka (Pvt.) Ltd. Tel: (94) 11 2785701 Fax: (94) 11 4410083-4	United Kingdom & Ireland 3M Health Care Ltd. Tel: 44-1509-611-611 Fax: 44-1509-237-288			
Chile 3M Chile S.A. Tel: (56-2) 4103000 Fax: (56-2) 4103400	Hungary 3M Hungária Kft. Tel: (36-1) 270-7777 Fax: (36-1) 320-0951	Middle East 3M Gulf Ltd. Tel: 971-4-3670-777 Fax: 971-4-3670-699	Spain 3M España, S.A. Tel: 34-91-321-6000 Fax: 34-91-321-6002	Venezuela 3M Manufacturera Venezuela Tel: 58-241-8391911 Fax: 58-241-8391928			
China 3M China Ltd. Tel: 86-21-6275-3535 Fax: 86-21-5208-2205	India 3M India Ltd. Tel: +91-80-2223 1414 Fax: +91-80-22231431	The Netherlands 3M Nederland B.V. Tel: 31-715-450-450 Fax: 31-715-450-212	Sweden 3M Svenska AB Tel: +46-8-92 21 00 Fax: +46-8-92 24 22	Viet Nam 3M Viet Nam Ltd. Tel: 84 8 640986 Fax: 84 8 641024			
Colombia 3M Colombia S.A. Tel: 57-1 4101666 Fax: 57-1 4161677	Indonesia PT 3M Indonesia Tel: (62-21) 520 3401 Fax: (62-21) 520 3106	New Zealand 3M New Zealand Tel: 64-9-477 4040 Fax: 64-9-477 6692	Switzerland & Eastern Europ 3M (Schweiz) AG Tel: +41 44 724 92 44 Fax: +41 44 724 94 80	e ISO 9001			

3M Peru S.A.

Philippines


3M Philippines

Tel: 511-224-2728

Fax: 511-224-3171

Tel: (632) 8133781 to 95

Fax: (632) 8145873 to 74



Tel.: 800 228-3957