

Technical Information Sheet Product Number 1776 3MTM Medical Nonwoven Tape

Effective: August 2013

Supersedes all previous specifications

Features and Benefits:

- Hypoallergenic Pressure Sensitive Adhesive - ETO and Gamma Sterilization Compatible

- Easy to Apply- Breathable- Conformable

Composition:

Tape caliper, without liner	~11.5 mils (0.3 mm) White Spunlace Polyester Nonwoven Tape		
Adhesive	Acrylate adhesive developed for medical/surgical use		
Backing	100% Polyester White Spunlace Nonwoven		
Liner	83 lb Poly-coated Kraft paper, silicone release on one side,		
	6 mils (0.15 mm)		

Typical Properties⁺:

Adhesion to Stainless Steel, 180 deg	25 ounces/inch width (0.7 kg/25mm width) (6.9 N/25 mm)
Liner Release	23 gm/inch width (23 gm/25 mm width)

⁺³M test methods are available upon request.

Roll Description: Tape supplied on 3 inch (76 mm) diameter cores.

Length	300 yards (274 m) ± 2 %
Width	48 inch (122 cm) <u>+</u> 1/32 in (1.0 mm)

Splices: Splices on the tape are intended for continuous processing.

Roll Width	Maximum Number of Splices per Roll
≤ 8 inch (20 cm)	3 per 100 yd (91 m)
	5 per 200 yd (183 m)
	9 per 300 yd (274 m)

Information on the splice frequency in roll widths greater than 8 inch (20 cm) is available upon request.

Packaging and Recommended Storage:

Product as supplied in original packaging will maintain stated properties for a period of two years from date stamped on shipping container when stored at temperatures of 50-80°F (10-27°C) and a relative humidity between 40-60 percent.

Please see reverse side for important product, safety and warranty information.

Product and Safety Information: User is solely responsible for determining the suitability of 3M samples and products for the intended use including any necessary safety or toxicity assessment. 3M will provide Material Safety Data Sheets and results of toxicity testing upon request. In every case before using any product in full scale production users should conduct their own tests to determine to their own satisfaction whether the product is of acceptable quality and is suitable for their particular purposes under their own operating conditions.

Notice: Nothing contained herein shall be construed to imply the nonexistence of any relevant patents or to constitute a permission, inducement or recommendations to practice any invention covered by any patent, without authority from the owners of this patent.

All statements, technical information and recommendations herein are based on tests 3M believes to be reliable, but the accuracy or completeness thereof is not guaranteed. 3M warrants only that products will meet 3M's specifications at the time of shipment to the customer. 3M does not offer any other warranty and does not warrant the performance, safety or such other characteristics of Products in combination with other materials. 3M specifically DOES NOT warrant Products for any intended or unintended uses (whether or not foreseeable); for compatibility or suitability with other components or compatibility with any methods of manufacture or conversion. The foregoing warranty is made in lieu of all other warranties, expressed or implied, including the implied warranties of merchantability, fitness for a particular purpose and freedom from non-infringement.

Limitation of Remedies: If products are proven not to meet 3M's specifications, the sole and exclusive remedy available and 3M's only obligation shall be, at 3M's option, to replace such quantity of Products which are proven out of specification or to refund the purchase price paid for Products.

Limitations of Liabilities: The remedies provided herein are exclusive remedies against 3M for any alleged or actual nonconformance to specifications or defect or other failure in products or for 3M's performance of its supply obligations. Under no circumstances is 3M liable for any direct, indirect, incidental, special or consequential damages (including lost profits) in any way related to the product under any theory of law including, but not limited to, negligence and strict liability.

Ordering Information

Call us at Customer Service to place an order: 800-742-1994 (U.S.). Visit our website: www.3M.com/medicalspecialties for product and services information, news about conferences we will be attending, new product highlights or to make a direct inquiry.

To have a sales representative contact you, to request samples or clinical and safety summaries, please contact us at 3M HELPLINE 800-228-3957 (U.S.) or for international inquiries, please contact your local country representative. Our 3M Medical Specialties subsidiary contacts are listed below.

International Locations							
Argentina 3M Argentina S.A.C.I.F.I.A. Tel.: 54-11-4339-2400 Fax: 54-11-4339-2640	Dominican Republic 3M Dominicana Tel: (809) 530-6560 X499 Fax: (809) 530-2960	Japan 3M Health Care Ltd. Tel: 81-3-3709-8289 Fax: 81-3-3709-8754	Poland 3M Poland Sp.z o.o. Tel: (48 22) 739 60 00 Fax: (48 22) 739 60 01	Thailand 3M Thailand Ltd. Tel: (66)2 260 8577 Ext.175 Fax: (66)2 261 7535			
Australia 3M Health Care Tel: 1 300 363 878 Fax: 1 800 060 888	Finland Suomen 3M Oy Tel.: 358-9-52-521 Fax: 358-9-512-29-44	Kazakhstan 3M Representation Office Tel: +7 3272 509 944 Fax: +7 3272 509 573	Portugal 3M Portugal Tel: 351-21-3134500 Fax: 351-21-3134680	Taiwan 3M Taiwan Ltd. Tel: 886-2-2704-9011 Fax: 886-2-2706-0355			
Belgium	France	Korea	Russia	Turkey			
3M Belgium N.V./S.A.	Laboratoires 3M Santé	3M Korea Ltd.	3M Russia	3M Turkey			
Tel: 32-2-722-5111	Tel: 33-1-3031-8376	Tel: 82-2-3771-4286	Tel: 7 (495) 784-74-74	Tel: 9 (0212) 3507777			
Fax: 32-2-720-0225	Fax: 33-1-3031-8378	Fax: 82-2-786-2825	Fax: 7 (495) 784-74-75	Fax: 9 (0212) 2821742			
Brazil 3M do Brasil Ltda. Tel: 55-19-3838-6661 Fax: 55-19-3838-7449	Germany	Malaysia	Singapore	Ukraine			
	3M Medica	3M Malaysia Sdn. Bhd.	3M Technologies (S) Pte. Ltd.	3M Ukraine			
	Tel: 49 2131 14 40 00	Tel: 6-03-7806 2888	Tel: 65-6450-8866	Tel: +380 44 490 57 77			
	Fax: 49 2131 14 49 99	Fax: 6-03-7806 2902	Fax: 65-6455-2130	Fax: +380 44 490 57 75			
Canada	Hong Kong	Mexico	Sri Lanka	United Kingdom & Ireland			
3M Canada	3M Hong Kong Ltd.	3M México, S.A. de C.V.	3M Lanka (Pvt.) Ltd.	3M Health Care Ltd.			
Tel: 800 364-3577	Tel: (852) 2806 6111	Tel: +52-55-5270-0400	Tel: (94) 11 2785701	Tel: 44-1509-611-611			
Fax: 800 341-4630	Fax: (852) 2234 6044	Fax: +52-55-5270-0433	Fax: (94) 11 4410083-4	Fax: 44-1509-237-288			
Chile	Hungary	Middle East	Spain	Venezuela			
3M Chile S.A.	3M Hungária Kft.	3M Gulf Ltd.	3M España, S.A.	3M Manufacturera Venezuela			
Tel: (56-2) 4103000	Tel: (36-1) 270-7777	Tel: 971-4-3670-777	Tel: 34-91-321-6000	Tel: 58-241-8391911			
Fax: (56-2) 4103400	Fax: (36-1) 320-0951	Fax: 971-4-3670-699	Fax: 34-91-321-6002	Fax: 58-241-8391928			
China 3M China Ltd. Tel: 86-21-6275-3535	India	The Netherlands	Sweden	Viet Nam			
	3M India Ltd.	3M Nederland B.V.	3M Svenska AB	3M Viet Nam Ltd.			
	Tel: +91-80-2223 1414	Tel: 31-715-450-450	Tel: +46-8-92 21 00	Tel: 84 8 640986			

Switzerland & Eastern Europe

3M (Schweiz) AG Tel: +41 44 724 92 44 Fax: +41 44 724 94 80

Fax: +46-8-92 24 22

Fax: 86-21-5208-2205 Colombia

3M Colombia S.A. Tel: 57-1 4101666 Fax: 57-1 4161677

Czech Republic

3M Cesko s.r.o. Tel: +420-261 380 111 Fax: +420-261 380 110

Denmark

3M a/sTel: 45-43-48-0100 Fax: 45-43-96-8596 Fax: +91-80-22231431

Indonesia

PT 3M Indonesia Tel: (62-21) 520 3401 Fax: (62-21) 520 3106

3M Israel Ltd. Tel: 972-9-9615000 Fax: 972-9-9615050

3M Italia SPA Tel: +39 2 7035-2406 Fax: +39 2 7035-2484 Fax: 31-715-450-212

New Zealand

3M New Zealand Tel: 64-9-477 4040 Fax: 64-9-477 6692

3M Peru S.A. Tel: 511-224-2728 Fax: 511-224-3171

Philippines

3M Philippines Tel: (632) 8133781 to 95 Fax: (632) 8145873 to 74

3M Medical Specialties

3M Center, Building 275-5W-05 St. Paul. MN 55144-1000 USA Tel.: 800 228-3957



Fax: 84 8 641024