

Technical Information Sheet Product Number 9874 3M™ Double Coated Medical Tape

Effective: June 2009

Supersedes all previous specifications

Features and Benefits:

Hypoallergenic pressure sensitive adhesive

- EtO & Gamma compatible

- Clear film

- Medium weight release liner
- Easy to apply
- Adheres to a wide range of substrates

Composition:

Tape caliper w/o liner	4.9 mil (0.12 mm) Double Coated Polyethylene Tape
Carrier	3 mil (0.08 mm) Transparent Polyethylene Film
Adhesive	Acrylate designed for medical/surgical use
Release Liner	60 lb. Bleached Kraft paper, silicone release coated both sides (3.6 mil/0.09 mm)

Typical Properties⁺:

180° Peel Adhesion to Stainless Steel	19 oz/inch width (0.54 Kg/25 mm) (5.3 N/25 mm)
Liner Release	8 gm/inch width (8 gm/25 mm)

⁺³M test methods are available upon request.

Roll Description: Tape supplied on 3 inch (76 mm) diameter cores.

Length	700 yards (640 m) ± 2% Maximum
Width	48 inch (122 cm) \pm 1/32 in (1.0 mm)

Splices: Splices on the tape are intended for continuous processing.

Roll Width	Maximum Number of Splices per Roll
≤ 8 inch (20 cm)	2 per 100 yd (91 m)
	3 per 200 yd (183 m)
	8 per 700 yd (640 m)

Information on the splice frequency in roll widths greater than 8 inch (20 cm) is available upon request.

Packaging and Recommended Storage

Product as supplied in original packaging will maintain stated properties for a period of two years from date stamped on shipping container when stored at temperatures of 50-80°F (10-27°C) and a relative humidity of 40-60 percent.

Please see reverse side for important product, safety and warranty information.

Product and Safety Information: User is solely responsible for determining the suitability of 3M samples and products for the intended use including any necessary safety or toxicity assessment. 3M will provide Material Safety Data Sheets and results of toxicity testing upon request. In every case before using any product in full scale production users should conduct their own tests to determine to their own satisfaction whether the product is of acceptable quality and is suitable for their particular purposes under their own operating conditions.

Notice: Nothing contained herein shall be construed to imply the nonexistence of any relevant patents or to constitute a permission, inducement or recommendations to practice any invention covered by any patent, without authority from the owners of this patent.

All statements, technical information and recommendations herein are based on tests 3M believes to be reliable, but the accuracy or completeness thereof is not guaranteed. 3M warrants only that products will meet 3M's specifications at the time of shipment to the customer. 3M does not offer any other warranty and does not warrant the performance, safety or such other characteristics of Products in combination with other materials. 3M specifically DOES NOT warrant Products for any intended or unintended uses (whether or not foreseeable); for compatibility or suitability with other components or compatibility with any methods of manufacture or conversion. The foregoing warranty is made in lieu of all other warranties, expressed or implied, including the implied warranties of merchantability, fitness for a particular purpose and freedom from non-infringement.

Limitation of Remedies: If products are proven not to meet 3M's specifications, the sole and exclusive remedy available and 3M's only obligation shall be, at 3M's option, to replace such quantity of Products which are proven out of specification or to refund the purchase price paid for Products.

Limitations of Liabilities: The remedies provided herein are exclusive remedies against 3M for any alleged or actual nonconformance to specifications or defect or other failure in products or for 3M's performance of its supply obligations. Under no circumstances is 3M liable for any direct, indirect, incidental, special or consequential damages (including lost profits) in any way related to the product under any theory of law including, but not limited to, negligence and strict

Ordering Information

Call us at Customer Service to place an order: 800-742-1994 (U.S.). Visit our website: www.3M.com/medicalspecialties for product and services information, news about conferences we will be attending, new product highlights or to make a direct inquiry.

To have a sales representative contact you, to request samples or clinical and safety summaries, please contact us at 3M HELPLINE 800-228-3957 (U.S.) or for international inquiries, please contact your local country representative. Our 3M Medical Specialties subsidiary contacts are listed below.

International Locations Argentina **Dominican Republic** Kazakhstan **Poland Thailand** 3M Argentina S.A.C.I.F.I.A. 3M Poland Sp.z o.o. 3M Dominicana 3M Representation Office 3M Thailand Ltd. Tel.: 54-11-4339-2400 Tel: (809) 530-6560 X499 Tel: +7 3272 509 944 Tel: (48 22) 739 60 00 Tel: (66)2 260 8577 Ext.175 Fax: 54-11-4339-2640 Fax: (809) 530-2960 Fax: +7 3272 509 573 Fax: (48 22) 739 60 01 Fax: (66)2 261 7535 Australia Finland **Portugal** Taiwan Korea 3M Health Care Suomen 3M Oy 3M Korea Ltd. 3M Portugal 3M Taiwan Ltd. Tel: 1 300 363 878 Tel.: 358-9-52-521 Tel: 82-2-3771-4286 Tel: 351-21-3134500 Tel: 886-2-2704-9011 Fax: 1 800 060 888 Fax: 358-9-512-29-44 Fax: 82-2-786-2825 Fax: 351-21-3134680 Fax: 886-2-2706-0355 **Belgium** France Malaysia Russia Turkey 3M Belgium N.V./S.A. Laboratoires 3M Santé 3M Malaysia Sdn. Bhd. 3M Russia 3M Turkey Tel: 6-03-7806 2888 Tel: 9 (0212) 3507777 Tel: 32-2-722-5111 Tel: 33-1-3031-8376 Tel: 7 (495) 784-74-74 Fax: 32-2-720-0225 Fax: 33-1-3031-8378 Fax: 6-03-7806 2902 Fax: 7 (495) 784-74-75 Fax: 9 (0212) 2821742 Brazil Germany Singapore Ukraine 3M México, S.A. de C.V. 3M do Brasil Ltda. 3M Technologies (S) Pte. Ltd. 3M Ukraine 3M Medica Tel: 55-19-3838-6661 Tel: 49 2131 14 40 00 Tel: +52-55-5270-0400 Tel: 65-6450-8866 Tel: +380 44 490 57 77 Fax: 55-19-3838-7449 Fax: 49 2131 14 49 99 Fax: +52-55-5270-0433 Fax: 65-6455-2130 Fax: +380 44 490 57 75 Canada Hong Kong Middle East Sri Lanka United Kingdom & Ireland 3M Canada 3M Hong Kong Ltd. 3M Gulf Ltd. 3M Lanka (Pvt.) Ltd. 3M Health Care Ltd. Tel: 800 364-3577 Tel: (852) 2806 6111 Tel: 971-4-3670-777 Tel: (94) 11 2785701 Tel: 44-1509-611-611 Fax: (852) 2234 6044 Fax: 971-4-3670-699 Fax: 44-1509-237-288 Fax: 800 341-4630 Fax: (94) 11 4410083-4 Hungary The Netherlands Venezuela

Chile

3M Chile S.A. Tel: (56-2) 4103000 Fax: (56-2) 4103400

China

3M China Ltd. Tel: 86-21-6275-3535 Fax: 86-21-5208-2205

Colombia

3M Colombia S.A. Tel: 57-1 4101666 Fax: 57-1 4161677

Czech Republic

3M Cesko s.r.o. Tel: +420-261 380 111 Fax: +420-261 380 110

Denmark

3M a/sTel: 45-43-48-0100 Fax: 45-43-96-8596 Indonesia PT 3M Indonesia

3M India Ltd.

3M Hungária Kft.

India

Tel: (36-1) 270-7777

Fax: (36-1) 320-0951

Tel: +91-80-2223 1414

Fax: +91-80-22231431

Tel: (62-21) 520 3401 Fax: (62-21) 520 3106

3M Italia SPA Tel: +39 2 7035-2406

3M Health Care Ltd. Tel: 81-3-3709-8289 Fax: 81-3-3709-8754

Fax: +39 2 7035-2484

3M New Zealand Tel: 64-9-444 4760 Fax: 64-9-443 7885

New Zealand

3M Nederland B.V.

Tel: 31-715-450-450

Fax: 31-715-450-212

Pern 3M Peru S.A. Tel: 511-224-2728 Fax: 511-224-3171

Philippines 3M Philippines Tel: (632) 8133781 to 95 Fax: (632) 8145873 to 74

Spain 3M España, S.A. Tel: 34-91-321-6000 Fax: 34-91-321-6002

Sweden 3M Svenska AB Tel: +46-8-92 21 00 Fax: +46-8-92 24 22

Switzerland & Eastern Europe

3M (Schweiz) AG Tel: +41 44 724 92 44 Fax: +41 44 724 94 80



3M Medical Specialties 3M Center, Building 275-5W-05 St. Paul. MN 55144-1000 USA Tel.: 800 228-3957



3M Manufacturera Venezuela

Tel: 58-241-8391911

Fax: 58-241-8391928

3M Viet Nam Ltd.

Tel: 84 8 640986

Fax: 84 8 641024

Viet Nam