

# **Technical Information Sheet** Product Numbers 9943 & 9943J $3M^{TM}$ Thin Hydrocolloid Adhesive Tape

Effective: July 2009

Supersedes all previous specifications

### **Features and Benefits:**

Very Gentle Adhesive, Good for Chronic/Repeat Wear - Conformable with Low Stretch

Good Adhesion, Short and Long Term (up to 4 days)

Translucent Pale Gold Color

Cushioning, Protective

- Breathable

- Gamma Compatible

- Water-Resistant Backing

# **Composition:**

Tape Caliper	19 mils (0.5 mm) without liner		
Backing	1 mil (0.03 mm) Polyurethane		
Adhesive	18 mil (0.5 mm) Hydrocolloid developed for medical/surgical use		
Release Liner	4.5 mil (0.11 mm) 80 lb. bleached Kraft paper, 2-side silicone		
	coated, Not Printable		

**Typical Properties:** 3M test methods are available upon request.

Absorbency, gms/gm weight, 48 hrs	3 grams buffered solution/1 gram adhesive		
Liner Release	30-80 gms/inch width (30-80 gms/25.4mm width)		

**Roll Description:** Tape supplied on 3 inch (76 mm) diameter cores.

#9943J Jumbo Roll Length, Maximum	75 yds. (68.6 m)	
#9943J Jumbo Tape Width, Range*	9.25–9.75 inches (23.5–24.8cm)	
Target	9.5"(24.1cm)	
#9943J Jumbo Liner Width, Max.*	10.25 in. (26.0 cm)	
#9943 Slit Roll Length, Maximum	1 inch (25mm) width – 75yds. (68.6 m)	
my to said from Bengun, nammum	>1 inch (25mm) width – 150yds. (137.1 m)	
#9943 Slit Tape Width, Maximum*	9 inches (22.9 cm) $\pm 1/16$ inch	
#9943 Slit Liner Width, Maximum *	9.25 inches (23.5 cm)	
Roll Diameter	Approx. 10 inches at 75yds. (25.4 cm-68.6m)	

\*Note: #9943J jumbo rollstock has extended liner on both edges and length is limited to 75 yards (68.6m). All #9943 slit rollstock is sold with 1/8" (3.2 mm) extended liner (EL) on both edges and up to 150 yards long.

**Splices:** Splices on the tape and liner are intended for continuous processing.

	<b>Maximum Number of Splices per Roll</b>	
3 ]	per 75 yards (68.6 m); 6 per 150 yds. (13	7.1 m)

### **Packaging and Recommended Storage**

Product as supplied in original packaging will maintain stated properties for a period of 2 years from date stamped on shipping container when stored at temperatures of 50-80°F (10-27°C) and relative humidity of 40-60 percent.

Please see reverse side for important product, safety and warranty information.

Product and Safety Information: User is solely responsible for determining the suitability of 3M samples and products for the intended use including any necessary safety or toxicity assessment. 3M will provide Material Safety Data Sheets and results of toxicity testing upon request. In every case before using any product in full scale production users should conduct their own tests to determine to their own satisfaction whether the product is of acceptable quality and is suitable for their particular purposes under their own operating conditions.

Notice: Nothing contained herein shall be construed to imply the nonexistence of any relevant patents or to constitute a permission, inducement or recommendations to practice any invention covered by any patent, without authority from the owners of this patent.

All statements, technical information and recommendations herein are based on tests 3M believes to be reliable, but the accuracy or completeness thereof is not guaranteed. 3M warrants only that products will meet 3M's specifications at the time of shipment to the customer. 3M does not offer any other warranty and does not warrant the performance, safety or such other characteristics of Products in combination with other materials. 3M specifically DOES NOT warrant Products for any intended or unintended uses (whether or not foreseeable); for compatibility or suitability with other components or compatibility with any methods of manufacture or conversion. The foregoing warranty is made in lieu of all other warranties, expressed or implied, including the implied warranties of merchantability, fitness for a particular purpose and freedom from non-infringement.

Limitation of Remedies: If products are proven not to meet 3M's specifications, the sole and exclusive remedy available and 3M's only obligation shall be, at 3M's option, to replace such quantity of Products which are proven out of specification or to refund the purchase price paid for Products.

Limitations of Liabilities: The remedies provided herein are exclusive remedies against 3M for any alleged or actual nonconformance to specifications or defect or other failure in products or for 3M's performance of its supply obligations. Under no circumstances is 3M liable for any direct, indirect, incidental, special or consequential damages (including lost profits) in any way related to the product under any theory of law including, but not limited to, negligence and strict

## **Ordering Information**

Call us at Customer Service to place an order: 800-742-1994 (U.S.). Visit our website: www.3M.com/medicalspecialties for product and services information, news about conferences we will be attending, new product highlights or to make a direct inquiry.

To have a sales representative contact you, to request samples or clinical and safety summaries, please contact us at 3M HELPLINE 800-228-3957 (U.S.) or for international inquiries, please contact your local country representative. Our 3M Medical Specialties subsidiary contacts are listed below.

International Locations						
	Argentina 3M Argentina S.A.C.I.F.I.A. Tel.: 54-11-4339-2400 Fax: 54-11-4339-2640	<b>Dominican Republic</b> 3M Dominicana Tel: (809) 530-6560 X499 Fax: (809) 530-2960	<b>Kazakhstan</b> 3M Representation Office Tel: +7 3272 509 944 Fax: +7 3272 509 573	<b>Poland</b> 3M Poland Sp.z o.o. Tel: (48 22) 739 60 00 Fax: (48 22) 739 60 01	<b>Thailand</b> 3M Thailand Ltd. Tel: (66)2 260 8577 Ext.175 Fax: (66)2 261 7535	
	Australia 3M Health Care Tel: 1 300 363 878 Fax: 1 800 060 888	Finland Suomen 3M Oy Tel.: 358-9-52-521 Fax: 358-9-512-29-44	<b>Korea</b> 3M Korea Ltd. Tel: 82-2-3771-4286 Fax: 82-2-786-2825	Portugal 3M Portugal Tel: 351-21-3134500 Fax: 351-21-3134680	<b>Taiwan</b> 3M Taiwan Ltd. Tel: 886-2-2704-9011 Fax: 886-2-2706-0355	
	Belgium 3M Belgium N.V./S.A. Tel: 32-2-722-5111 Fax: 32-2-720-0225	France Laboratoires 3M Santé Tel: 33-1-3031-8376 Fax: 33-1-3031-8378	Malaysia 3M Malaysia Sdn. Bhd. Tel: 6-03-7806 2888 Fax: 6-03-7806 2902	<b>Russia</b> 3M Russia Tel: 7 (495) 784-74-74 Fax: 7 (495) 784-74-75	<b>Turkey</b> 3M Turkey Tel: 9 (0212) 3507777 Fax: 9 (0212) 2821742	
	Brazil 3M do Brasil Ltda. Tel: 55-19-3838-6661 Fax: 55-19-3838-7449	Germany 3M Medica Tel: 49 2131 14 40 00 Fax: 49 2131 14 49 99	<b>Mexico</b> 3M México, S.A. de C.V. Tel: +52-55-5270-0400 Fax: +52-55-5270-0433	<b>Singapore</b> 3M Technologies (S) Pte. Ltd. Tel: 65-6450-8866 Fax: 65-6455-2130	Ukraine 3M Ukraine Tel: +380 44 490 57 77 Fax: +380 44 490 57 75	
	<b>Canada</b> 3M Canada Tel: 800 364-3577 Fax: 800 341-4630	Hong Kong 3M Hong Kong Ltd. Tel: (852) 2806 6111 Fax: (852) 2234 6044	Middle East 3M Gulf Ltd. Tel: 971-4-3670-777 Fax: 971-4-3670-699	<b>Sri Lanka</b> 3M Lanka (Pvt.) Ltd. Tel: (94) 11 2785701 Fax: (94) 11 4410083-4	<b>United Kingdom &amp; Ireland</b> 3M Health Care Ltd. Tel: 44-1509-611-611 Fax: 44-1509-237-288	
	Chile	Hungary	The Netherlands	Spain	Venezuela	

Chile

3M Chile S.A. Tel: (56-2) 4103000 Fax: (56-2) 4103400

China

3M China Ltd. Tel: 86-21-6275-3535 Fax: 86-21-5208-2205

Colombia

3M Colombia S.A. Tel: 57-1 4101666 Fax: 57-1 4161677

Czech Republic

3M Cesko s.r.o. Tel: +420-261 380 111 Fax: +420-261 380 110

Denmark

3M a/sTel: 45-43-48-0100 Fax: 45-43-96-8596 Fax: +91-80-22231431 Indonesia

3M India Ltd.

India

3M Hungária Kft.

Tel: (36-1) 270-7777

Fax: (36-1) 320-0951

Tel: +91-80-2223 1414

PT 3M Indonesia Tel: (62-21) 520 3401 Fax: (62-21) 520 3106

3M Italia SPA Tel: +39 2 7035-2406

3M Health Care Ltd. Tel: 81-3-3709-8289 Fax: 81-3-3709-8754

Fax: +39 2 7035-2484

Tel: 64-9-477 4040 Fax: 64-9-477 6692

New Zealand

3M New Zealand

Pern 3M Peru S.A. Tel: 511-224-2728

Fax: 511-224-3171

3M Nederland B.V.

Tel: 31-715-450-450

Fax: 31-715-450-212

**Philippines** 3M Philippines

Tel: (632) 8133781 to 95 Fax: (632) 8145873 to 74 3M Manufacturera Venezuela Tel: 58-241-8391911 Fax: 58-241-8391928

Viet Nam 3M Viet Nam Ltd. Tel: 84 8 640986 Fax: 84 8 641024

Switzerland & Eastern Europe

3M (Schweiz) AG Tel: +41 44 724 92 44 Fax: +41 44 724 94 80

3M España, S.A.

3M Svenska AB

Tel: +46-8-92 21 00

Fax: +46-8-92 24 22

Sweden

Tel: 34-91-321-6000

Fax: 34-91-321-6002



**3M Medical Specialties** 3M Center, Building 275-5W-05 St. Paul. MN 55144-1000 USA

Tel.: 800 228-3957

